EN: NECKSTRECH – USER MANUAL

Gravity Finger Massager


According to the relevant data, there are 6.5 million people suffering from cervical spondylosis. Take care of your body and maintain correct posture. Restore your vertebrae to your natural position with ergonomic support, relax your muscles, relieve stress, or tension and relax.

Innovative gravity finger massage

Under the effect of gravity, the posture of the cervical spine (vertebrae) can be effectively improved, restoring your neck's natural curvature, eliminating nerve compression and insufficient blood flow caused by dislocation.

There are three types of physiological curvature manifestations of cervical vertebrae.

1. Type "C" (Normal/natural curvature)
2. Straightened (Lessened curvature) and
3. Reverse bow (Severely lessened curvature)

USE METHOD

Use this product if you feel sore or experiencing neck pain, tension, tingling, weak arms, numb fingers, bulging pain in your shoulders and back, dizziness and nausea.

1. Rest in a supine position.
2. Position NECKSTRETCH below your neck.
3. Bend your knees and extend your hands to the front or back.
4. Your head is positioned about 1 cm above the surface.

Neck Traction Massage

Place the product on a flat surface, lie on your back, relax your neck naturally, and allow gravity to stretch and massage your head.

Neck Muscle Massage

Fix the product in the right position, lie on your back, shake your neck to the left and right, and you can play the neck muscle.

Other methods

Hold the massager handle and place it on the thigh, shoulder and other curved protruding parts of the human body to knead and massage them.

Note: please refrain from using excessive force on your neck. Do not move your head with excessive frequency or speed.

Correct posture

1. Concave surface: Head suspension, strong downward pressure
2. Convex surface: The head is supported with a light pressure.

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